	DE	Gebrauchsanweisung LCU Hüftprothesensystem
	EN	Instructions for use LCU Hip System
	BG	Инструкции за употреба система за тазобедрена става LCU
	CS	Návod k použití Kýčelní systém LCU
	DA	Brugsanvisning LCU-hoftesystem
	EL	Οδηγίες χρήσης σύστημα ισχίου LCU
	ES	Instrucciones de uso Sistema de cadera LCU
	ET	Kasutusjuhend LCU puusasüsteem
	FI	Käyttöohje LCU-lonkkajärjestelmä
	FR	Mode d'emploi Système de prothèse de hanche LCU
	HR	Upute za korištenje Sustav kuka LCU
	HU	Használati útmutató LCU csípőrendszer
	ID	Instruksi Penggunaan Sistem Pinggul LCU
	IT	Istruzioni per l'uso Sistema per anca LCU
	LT	Naudojimo instrukcijos LCU klubo sistema
	LV	Lietošanas instrukcija LCU gūžas sistēma
	NL	Gebruiksaanwijzing LCU heupsysteem
	NO	Bruksanvisning LCU-hoftesystemet
	PL	Instrukcja użycia System biodra LCU
	PT	Instruções de utilização Sistema de anca LCU
	RO	Instructiuni de folosire Sistemul LCU pentru șold
	RU	Инструкция по применению Имплантаты Тазобедренный сустав LCU
	SK	Inštrukcie na používanie Systém bedrového kĺbu LCU
	SL	Navodila za uporabo Kolčni sistem LCU
	SV	Bruksanvisning Implantat LCU höftsystem
	TR	Kullanım Talimatı İmplantlar LCU Kalça Sistemi
	UK	Інструкція для застосування система для ендопротезування кульшового суглоба LCU
	ZH	说明书 骨科植入物 LCU 髋关节系

Not for US distribution.

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CN DIST

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China Distributor, After Sales Service

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Insertion Forceps				
134-141/00	-	-	-	M12
Resection Guide				
175-310/05	-	-	-	M13
Instrument Tray				
165-100/10	165-100/11	-	-	M8, M10, M36
Head Impactor, with exchangeable plastic head				
175-360	-	-	-	M10, M15, M21, M23

REF	LCU Hip System – Instrument Set			CE	MAT
X-Ray Template					
165-140/35	165-141/35	165-142/35	165-143/35	M24	
165-144/35	165-150/00	-	-		

REF	LCU Hip System – Instrument Set			CE ₀₄₈₂	MAT
Taper Cap					
179-122/01	-	-	-	M8	
Rasp Handle, with quick coupling					
130-394/01	-	-	-	M1, M12, M14, M15	
Handle, with quick coupling					
130-394/02	130-394/03	-	-	M1, M12, M14, M15	
Positioning Guide, for alignment of anteversion					
130-393/81	-	-	-	M15	
Bone Plug Impactor					
131-220	-	-	-	M11, M29	
Impactor					
130-622/01	-	-	-	M12, M13, M15, M37	
Instrument Set complete					
165-100/30	165-100/31	-	-	Refer to single article number	
Positioner					
130-711	-	-	-	M15	
Box Chisel					
130-716	-	-	-	M10, M13, M15	
Impactor					
130-713	-	-	-	M15	
Plastic Trial Head					
132-926/03	132-928/03	132-932/03	132-936/03	M23	
132-926/02	132-928/02	132-932/02	132-936/02		
132-928/04	132-932/04	132-936/04	132-926/01		
132-928/01	132-932/01	132-936/01	-		
Bone Compressor					
165-111/10	165-111/11	165-111/12	165-111/13	M17	
165-111/14	165-111/15	165-111/16	165-111/18		
165-111/20	165-111/08	165-111/09	-		
Trial Neck Segment					
165-110/25	165-110/30	-	-	M15	
Femoral Canal Opener					
130-617	-	-	-	M12, M15, M37	

MAT	Legende zu den Werkstoffen Λεγέδα на материалите Oversigt over materialer Leyenda de materiales Materiaalien selitykset Legenda materijala Legenda material Medžiagų išaiškinimas Legenda van materialen Legenda materialów Legenda materialelor Legenda k materiálov Förklaring av material Розшифрування матеріалів	Legend of materials Popisky k materiálům Υπόμνημα υλικών Materjalide loend Légende des matériaux Jelmagyarázat az anyagokhoz Legenda dei materiali Skaidrojums par materiāliem Forklaring for materialene Legenda dos materiais Расшифровка материалов Legenda materialov Materiyal Yazıları 材料图例
M1	Kobaltbasis-Legierung Сплав на основата на кобалт Kobolt-baseret legering Aleación a base de cobalto Kobolttipohjainen seos Legura na bazi kobalta Alói berbasis kobalt Kobalto pagrindo lydinys Kobaltlegering Stop na bazie kobaltu Aliaj pe bază de cobalt Zlitina na báze kobaltu Kobolt-baseret legering Сплав на основі кобальту CoCrMo, ISO 5832-12; ASTM F1537	Cobalt-based alloy Slitina na bázi kobaltu Κράμα με βάση κοβάλτιο Koobaltipõhine sulam Alliage à base de cobalt Kobaltalápú ötvözet Lega di cobalto Kobalta bāzes sakausējums Koboltbasert legering Liga à base de cobalto Сплав на основе кобальта Zlitina na osnovi kobalta Kobalt tabanlı alaşım 钴基合金
M2	Titan-Aluminium-Legierung Τιτανο-αλουμινεία σπлав Titan-aluminium-legering Aleación de titanio-aluminio Titaanin ja alumiinin seos Legura na bazi titanija Alói titanium aluminium Titano ir aluminio lydinys Titanium-aluminiumlegering Stop tytanu i aluminium Aliaj de titan și aluminii Zlitina titānu a hlinīku Titan-aluminiumlegering Τιταν-αλουμινέβιυ σπлав Ti6Al4V, DIN EN ISO 5832-3; ASTM F136	Titanium-aluminum alloy Slitina titanu a hlinīku Κράμα τιτανίου-αλουμινίου Titaani-/alumiiniumisulam Alliage titane-aluminium Titán-aluminium ötvözet Lega di titanio-alluminio Titāna-alumīnija sakausējums Titan-aluminium-legering Liga de titānio/alumīnio Τιτανο-αλουμινέβιυ σπлав Zlitina titana in aluminija Titanium-aluminiumyum alaşım 钛铝合金
M6	Calcium-Phosphat-Beschichtung Ποκρίτιε от калциев фосфат Calciumfosfatbelægning Recubrimiento de fosfato de calcio Kalsiumfosfaatpinnoite Premaz od kalcijevog fosfata Salut fosfat kalsium Kalcio fosfato danga Coating van calciumfosfaat Powłoka z fosforanu wapnia Ínvelis din fosfat de calciu Povrchová úprava s fosforečnanom vápenatým Kalciumfosfatbelægning Ποκρίτιττз фосфату кальцію ASTM F-1609, (HX)	Calcium phosphate coating Povlak z fosforečnanu vápenatého Επκόλυση φωσφορικού ασβεστίου KalsiumfosfaatKate Revêtement de phosphate de calcium Kalcium-foszfát bevonat Rivestimento di fosfato di calcio Kalcija fosfāta pārklājums Belegg av kalsiumfosfat Rivestimento de fosfato de cálcio Кальций-фосфатное покрытие Prevlēka iz kalcijevoga fosfata Kalsiyum fosfat kaplama 磷酸钙涂层
M8	Polyphenylsulfon Πολυφенилсульфон Polyfenylsulfon Polifenilsulfona Polyfenylisulfoni Polifenilsulfon Polifenylosulfon Polifenilsulfonas Polyfenylsulfon Polifenylosulfon Polifenilsulfonā Polyfenylsulfón Polyfenylsulfon Πολυφенилсульфон PPSU	Polyphenylsulfone Polyfenylsulfon Πολυφαινυλοσουλφόνη Polüfenülsulfoon Polyphénylsulfone Polifenilszulfon Polifenilsulfone Polifenilsulfons Polyfenylsulfon Polifenilsulfona Πολυφенилсульфон Polifenilsulfon Polifenilsulfon 聚苯砜

M10	Edelstahl Неръждаема стомана Rustfritt stål Acero inoxidable Ruostumaton teräs Nehrdajúci čelik Baja tahan karat Nerūdijantis plienas Roestvrij staal Stal nierdzewna Ořel inoxidable Nehrdzavejúca oceľ Rostfritt stål Нержавиоуча сталь X5CrNi18-10, DIN EN 10088-1	Stainless steel Nerezová ocel Ανοξειδωτός χάλυβας Roostevaba teras Acier inoxydable Rozsdamentes acél Acciaio inossidabile Nerūsējošais tērauds Rustfritt stål Aço inoxidável Нержавеющая сталь Nerjavno jeklo Paslanmaz çelik 不锈钢
M11	Edelstahl Неръждаема стомана Rustfritt stål Acero inoxidable Ruostumaton teräs Nehrdajúci čelik Baja tahan karat Nerūdijantis plienas Roestvrij staal Stal nierdzewna Ořel inoxidable Nehrdzavejúca oceľ Rostfritt stål Нержавиоуча сталь X8CrNiS 18-9, DIN EN 10088-1	Stainless steel Nerezová ocel Ανοξειδωτός χάλυβας Roostevaba teras Acier inoxydable Rozsdamentes acél Acciaio inossidabile Nerūsējošais tērauds Rustfritt stål Aço inoxidável Нержавеющая сталь Nerjavno jeklo Paslanmaz çelik 不锈钢
M12	Edelstahl Неръждаема стомана Rustfritt stål Acero inoxidable Ruostumaton teräs Nehrdajúci čelik Baja tahan karat Nerūdijantis plienas Roestvrij staal Stal nierdzewna Ořel inoxidable Nehrdzavejúca oceľ Rostfritt stål Нержавиоуча сталь X20Cr13, DIN EN 10088-1	Stainless steel Nerezová ocel Ανοξειδωτός χάλυβας Roostevaba teras Acier inoxydable Rozsdamentes acél Acciaio inossidabile Nerūsējošais tērauds Rustfritt stål Aço inoxidável Нержавеющая сталь Nerjavno jeklo Paslanmaz çelik 不锈钢
M13	Edelstahl Неръждаема стомана Rustfritt stål Acero inoxidable Ruostumaton teräs Nehrdajúci čelik Baja tahan karat Nerūdijantis plienas Roestvrij staal Stal nierdzewna Ořel inoxidable Nehrdzavejúca oceľ Rostfritt stål Нержавиоуча сталь X10CrNi 18-8, DIN EN 10270-3	Stainless steel Nerezová ocel Ανοξειδωτός χάλυβας Roostevaba teras Acier inoxydable Rozsdamentes acél Acciaio inossidabile Nerūsējošais tērauds Rustfritt stål Aço inoxidável Нержавеющая сталь Nerjavno jeklo Paslanmaz çelik 不锈钢
M14	Edelstahl Неръждаема стомана Rustfritt stål Acero inoxidable Ruostumaton teräs Nehrdajúci čelik Baja tahan karat Nerūdijantis plienas Roestvrij staal Stal nierdzewna Ořel inoxidable Nehrdzavejúca oceľ Rostfritt stål Нержавиоуча сталь X90CrMoV18, DIN EN 10088-3	Stainless steel Nerezová ocel Ανοξειδωτός χάλυβας Roostevaba teras Acier inoxydable Rozsdamentes acél Acciaio inossidabile Nerūsējošais tērauds Rustfritt stål Aço inoxidável Нержавеющая сталь Nerjavno jeklo Paslanmaz çelik 不锈钢

M15	Edelstahl Нерждаема стомана Rustfritt stål Acero inoxidable Ruostumatón teräs Nehrdajuci čelik Baja tahan karat Nerūdijantis plienas Roestvrij staal Stal nierdzewna Ođel inoxidabil Nehrdzavejūca ocel' Rostfritt stål Нержавіюча сталь X5CrNiCuNb 16-4 / X5CrNiCuNb 17-4, ISO 7153-1, AISI 630	Stainless steel Nerezová ocel Ανοξειδωτός χάλυβας Roostevaba teras Acier inoxydable Rozsdamentes acél Acciaio inossidabile Nerūsējošais tērauds Rustfritt stål Aço inoxidável Нержавеющая сталь Nerjavno jeklo Paslanmaz çelik 不锈钢
M17	Edelstahl Нерждаема стомана Rustfritt stål Acero inoxidable Ruostumatón teräs Nehrdajuci čelik Baja tahan karat Nerūdijantis plienas Roestvrij staal Stal nierdzewna Ođel inoxidabil Nehrdzavejūca ocel' Rostfritt stål Нержавіюча сталь X17CrNi16-2, DIN EN 10088-3, AISI 431	Stainless steel Nerezová ocel Ανοξειδωτός χάλυβας Roostevaba teras Acier inoxydable Rozsdamentes acél Acciaio inossidabile Nerūsējošais tērauds Rustfritt stål Aço inoxidável Нержавеющая сталь Nerjavno jeklo Paslanmaz çelik 不锈钢
M21	Silikon Силиций Silicone Silicona Silikoni Silikon Silikon Silikonas Siliconen Silikon Silicon Silikon Silikon Силикон BGA XV	Sillicone Silikon Σιλίκόνη Silikoon Silicone Szikon Silicone Silikons Silikon Silicone Силикон Silikon Silikon 硅酮
M23	Polypropylen-Номopolimer Полипропилен хомополимер Polypropylen-homopolymer Номopolimero de polipropileno Polypropyleenihomopolymeeri Номopolimer polipropilena Polipropilena homopolimer Polipropileno homopolimeras Polypropyleen homopolymeer Номopolimer polipropylenu Polipropilenā homopolimer Polypropylénový homopolymér Polypropren-homopolymer Полипропилену гомополимер PP-H, DIN EN 15860, ASTM D4101	Polypropylene homopolymer Polypropylen homopolymer Ομοπολυμερές πολυπροπυλένιο Polüpropüleenhomopolümeer Polypropylène homopolymère Polipropilén homopolimer Polipropilene omopolimero Polipropilēna homopolimērs Polypropylen-homopolymer Номopolimero de polipropileno Гомополимер полипропилена Номopolimer polipropilena Polipropilen homopolimeri 聚丙烯均聚物
M24	Polyvinylchlorid Поливинил хлорид Polyvinyl-klorid Cloruro de polivinilo Polyvinylikloridi Poli(vinil-klorid) Polivinil klorida Polivinilchloridas Polyvinylchloride Polichlorek winylu Clorurā de polivinil Polyvinylchlorid Polyvinyliklorid Полівінілхлорид PVC	Polyvinyl chloride Polyvinylchlorid Πολυβινυλοχλωρίδιο Polüvinüükloriid Chlorure de polyvinyle Polivinilklorid Cloruro di polivinile Polivinilchlorids Polyvinylklorid Policloreto de vinil Поливинилхлорид Polivinil klorid Polivinil klorür 聚氯乙烯

M29	Edelstahl Нерждаема стомана Rustfritt stål Acero inoxidable Ruostumatón teräs Nehrdajuci čelik Baja tahan karat Nerūdijantis plienas Roestvrij staal Stal nierdzewna Ođel inoxidabil Nehrdzavejūca ocel' Rostfritt stål Нержавіюча сталь X2CrNiMo 18-15-3, DIN ISO 5832-1:97, ASTM F138-00, ASTM F139-00	Stainless steel Nerezová ocel Ανοξειδωτός χάλυβας Roostevaba teras Acier inoxydable Rozsdamentes acél Acciaio inossidabile Nerūsējošais tērauds Rustfritt stål Aço inoxidável Нержавеющая сталь Nerjavno jeklo Paslanmaz çelik 不锈钢
M36	Thermoplastisches Elastomer Термопластичен еластомер Termoplastisk elastomer Elastómero termoplástico Termoplastinen elastomeeri Termoplastični elastomer Elastomer termoplastik Termoplastinis elastomeras Thermoplastisch elastomeer Elastomer termoplastyczny Elastomer termoplastic Termoplastický elastomér Termoplastiska elastomerer Термопластичний еластомер TPE	Thermoplastic elastomer Termoplastický elastomer Θερμοπλαστικό ελαστομερές Termoplastne elastomeer Chlorure de polyvinyle Hőre lágyuló elasztomer Elastomero termoplastico Termoplastisks elastomērs Termoplastelastomer Elastómero termoplástico Термопластичный эластомер Termoplastični elastomer Termoplastik elastomer 热塑性弹性体
M37	Silikonkautschuk Силициев каучук Silikonegummi Goma de silicona Silikonikumi Silikonska guma Karet silikon Silikoninė guma Siliconenrubber Guma silikonowa Cauciuc siliconic Silikonová guma Silikongummi Силиконовый каучук FDA§ 177.2600	Silicone Rubber Silikonová guma Ελαστικό σιλικόνης Silikoonkummi Caoutchouc silicone Szikilongumi Gomma siliconica Silikona gumija Silokongummi Borracha de silicone Силиконовая резина Silikonka guma Silikon Kauçuk 硅橡胶
M50	Kobaltbasis-Legierung Сплав на основата на кобалт Kobolt-baseret legering Aleación a base de cobalto Kobolttipohjainen seos Legura na bazi kobalta Aloi berbasis kobalt Kobalto pagrindo lydiny Kobaltlegering Stop na bazie kobaltu Aliaj pe bazā de cobalt Zlitina na báze kobaltu Kobolt-baseret legering Сплав на основі кобальту CoCrMo, ISO 5832-12; ASTM F799	Cobalt-based alloy Slitina na bázi kobaltu Κράμα με βάση κοβάλτιο Koobaltipõhine sulam Alliage à base de cobalt Kobaltalapú ötvözet Lega di cobalto Kobalta bāzes sakausējums Kobaltbasert legering Liga à base de cobalto Сплав на основе кобальта Zlitina na osnovi kobalta Kobalt tabanlı alaşım

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1. General

Please read this document carefully before using the system and keep it for future consultation! Failure to comply with these instructions for use may cause damage to patient's health, implants and instruments.

This document does not contain all of the information necessary for the selection and application of the system. For safe and proper handling, refer to further product-related instructions, such as the surgical technique associated with the system, the packaging and cleaning instructions and the device labels on the packaging. For instruments that require calibration, please refer to the Additional Product Information, Torque Wrench (IFU-EU-134-140-001).

2. User Group and Environment

The products may only be used and operated in an aseptic medical environment by persons who have the appropriate training, knowledge or experience in the orthopaedic and surgical field.

The intended users are physicians, theatre nurses as well as CSSD staff (CSSD = Central Sterile Supply Department).

3. Patient Group

The patient group for our medical device is made up of adult, anesthetized patients of any ethnic origin and of any gender, in whom one or more of the described indications are present and who are not considered unsuitable on the grounds of the listed contraindications.

4. System Description

The LCU Hip System is part of a mechanical reconstruction of the hip joint.

The LCU forms a hemiarthroplasty of the hip joint when combined with trauma heads (e.g. Modular Trauma Heads or Large Heads) and forms a total replacement of the hip joint if combined with prosthesis heads and an acetabular cup.

There are two necessary Instrument Sets to Implant LCU and the LCU cemented Hip System and some additional Instruments.

Instruments for Direct Anterior Approach: The instruments are designed for usage within the application of the Direct Anterior Approach.

5. Intended Use

Implants

The non-active, surgically-invasive implantable LCU Hip System manufactured by Waldemar Link GmbH & Co. KG is intended for long-term replacement of the femoral side of a diseased and / or defective hip joint in the human body. The LCU Hip System forms a replacement of the hip joint when combined with the prosthesis head and if applicable the acetabular cup. The LCU Hip System can be used in full-grown, anesthetized patients of any ethnic origin and sex. The LCU Hip System is implanted without cement (LCU HX, LCU PoroLink) and with cement (LCU Cemented).

The implants may only be used and operated in an aseptic medical environment by persons who have the required training, knowledge and experience in the orthopedic and surgical field. The implants are supplied in sterile condition individually packed as single-use products.

Instruments LCU Hip System

The purpose of the instrumentation is to allow the user to use the associated implant system as part of the procedures described in the associated surgical technique. Any other use of the instrumentation is not permitted.

The instruments consist of defined, combinable instruments.

All instruments of the instrument set are intended for transient use.

Instruments for Direct Anterior Approach

The purpose of the instrumentation is to allow the user to use the associated implant system as part of the procedures described in the associated surgical technique. Any other use of the instrumentation is not permitted.

The instruments consist of defined, combinable instruments.

All instruments of the instrument set are intended for transient use.

6. Indications

Implants

General Indications:

- Mobility-limiting diseases, fractures or defects of the hip joint or proximal femur which cannot be treated by conservative or osteosynthetic procedures

Indications:

- Primary and secondary osteoarthritis
- Rheumatoid arthritis
- Correction of functional deformities
- Avascular necrosis
- Femoral neck fractures
- Revision after implant loosening dependent on bone mass and quality

Instruments LCU Hip System

- The LCU and LCU cemented Hip System Instrument Set in the present combination is exclusively intended for the application and implantation of LCU and LCU cemented Hip System.

Instruments for Direct Anterior Approach

- The Direct Anterior Approach instrument set arrangement in the present combination is intended for the application and implantation of various hip reconstruction prostheses from LINK with the Direct Anterior Approach.

7. Contraindications

Implants

- Acute and chronic infections, local and systemic, insofar as they may compromise the successful implantation of a hip prosthesis (preoperative microbiological analysis recommended)
- Allergies to (implant) materials
- Insufficient / inadequate bone mass- or quality which prevents a stable anchorage of the prosthesis

Instruments LCU Hip System

- The LCU and LCU cemented Hip System Instrument Set in the present combination is not intended for the application and implantation of another Medical Device aside from LCU and LCU cemented Hip System.
- Material intolerance to the instrument materials

Instruments are not allowed to be implanted.

Instruments for Direct Anterior Approach

- The Direct Anterior Approach instrument set arrangement in the present combination is not intended for the application and implantation of another Medical Device aside from hip reconstruction prostheses from the LINK group.
- Material intolerance to the instrument materials

Instruments are not allowed to be implanted.

8. Possible Risks and Side Effects

Cemented

- Implant damage, implant fracture
- Infection
- Instability, Dislocation
- Malalignment
- Periprosthetic fracture
- Residual complaints
- Septic, aseptic loosening
- Subsidence
- Wear

Uncemented

- Aseptic Loosening
- Bone union problems
- Implant damage, implant fracture
- Infection
- Instability, Dislocation
- Malalignment
- Noises (e.g. cracking, popping, snapping, squeaking or grinding)
- Periprosthetic fracture
- Prosthetic impingement
- Residual complaints
- Soft tissue problems
- Stress Shielding
- Subsidence
- Wear

9. Clinical Benefit

The clinical benefit of the LCU Hip System is defined as:

- Relieve of pain, compared to untreated diseased hip joint
- Increased joint functionality, compared to untreated diseased hip joint
- Increased patient mobility, compared to untreated diseased hip joint

The clinical benefit of the instruments for the LCU and LCU cemented Hip System is the implantation of the LCU and LCU cemented Hip System.

The clinical benefit of the Direct Anterior Approach Instruments is the implantation of various hip reconstruction prostheses from LINK via the Direct Anterior Approach.

10. Implant Materials

Please refer to the relevant surgical technique associated with the system and the identification on the packaging for further information on implant materials. Further information on the material compositions is available from the manufacturer upon request.

- Cobalt-based alloy, CoCrMo, ISO 5832-12; ASTM F799
- Titanium-aluminum alloy (Ti6Al4V) according to ISO 5832-3 / ASTM F-136
- Calcium phosphate coating (HX) according to ASTM F-1609

11. Instrument Materials

Please refer to the REF list in this document and to the relevant surgical technique associated with the system and the identification on the packaging for further information on instrument materials.

Further information on the material compositions is available from the manufacturer upon request.

12. CMR Substances

Some system components contain cobalt as an alloy ingredient in a concentration above 0.1 % weight by weight.

Cobalt is listed as a CMR (carcinogenic, mutagenic and toxic to reproduction) substance.

The hazard class and category code(s) for cobalt are:

- Carc. 1B
- Repro. 1B

For identification of the affected components and further material specification, please refer to the REF list and the legend of materials in this document.

13. Implant Selection, permissible Combinations

For definitive identification information on the product such as system compatibility, article number, material and shelf life, refer to the identification on the implant and / or the packaging.

Also, refer to the following sections of this document and to the relevant surgical technique associated with the system for further information on implant selection and permissible combination and for information on allocating and handling the instruments to be used for the implantation.

Combinations with implants from other manufacturers and / or combinations with LINK implants that deviate from the surgical technique specifications have not been tested and are not permitted.

14. Permissible materials for tribological pairings of implants

The implant components of the system only form tribological pairings when combined with other Link implants, e. g. prosthesis heads.

When combined with other Link implants, the following tribological pairings may be formed:

- Titanium-aluminum alloy / CoCrMo alloy
- Titanium-aluminum alloy / Titanium-aluminum alloy
- Titanium-aluminum alloy / Composite ceramics (e.g. Aluminium oxide)

“Ceramic on metal” and “metal on metal” pairings are prohibited within the area of the bearing surface. The same applies for the pairing of BIOLOX delta with a ceramic component by another manufacturer.

15. Implant Size

Please refer to the relevant surgical technique associated with the system for further information on implant sizes.

16. Implant Anchoring

Implant components are labelled as to whether they are to be cemented or not.

17. Lifetime

The lifetime of our implants is limited in principle and is determined by individual factors such as, for example, body weight and the level of activity of the patient, as well as by the quality and professional execution of the implantation. Based on these individual influencing factors, Waldemar Link defines the overall average lifetime of an implant based on its survival rate (i.e. the proportion of functional implants after a certain period of time starting from the time of implantation). According to the results of the tests performed, the survival rate of our implants corresponds to the general state of the art at the time of approval of the implants.

The expected life time of instruments depends on material, design, application and processing. The expected life time of instruments by Waldemar Link GmbH & Co. KG is limited by restrictions in their usability and / or functionality.

18. Reprocessing / Reuse

The implants are supplied as sterile single-use devices. Implants whose protective packaging is opened or damaged, or implants which have already been implanted are not permitted to be reprocessed or reused.

The following risks may occur if implants are re-used:

- Infections
- Reduced implant lifetime
- Increased wear and wear debris complications
- Disease transmission
- Inadequate implant fixation
- Limited implant function
- Implant response and / or rejection

Instruments must be disinfected and sterilized prior to use. For more information please refer to the related chapters in this document and to the description in the reprocessing instructions H50.

Additionally observe our separate packaging and cleaning instructions for instruments.

Single-use products may not be reused.

19. Resterilization

The Implants are designed for single-use only. Resterilization is not permitted. Implants, as well as their materials are not suitable to be resterilized.

Unpredictable degradations may occur in these implants during resterilization. For sterilization of the Instruments, please refer to the description in the reprocessing instructions H50.

20. Storage and Transportation

Sterile-packaged implants must be stored in the undamaged original packaging in buildings with adequate protection against damage due to impacts, frost, humidity, excessive heat, and direct sunlight.

For Storage and Transport of the Instruments, please refer to the description in the reprocessing instructions H50.

21. Information for Patient Advisory

If the implantation of this system is considered to be the best solution for the patient and one of the circumstances described in section 22 is applicable to the patient, it is necessary to advise the patient with regard to the anticipated effects that these circumstances could have on the success of the operation. It is further recommended that the patient be informed about measures that he or she can take to reduce the effects of such complications. All information provided to the patient should be documented in writing by the operating surgeon. An implant ID for the patient must be handed over by the surgeon or hospital and the patient must be informed about the availability of a special patient information.

The patients should also be instructed:

- in detail about the surgery-related risks.
- in detail about the limitations of the implants, especially about the effects of excessive stress caused by body weight and physical activity, among other things. They should be encouraged to adjust their activities accordingly.
- about possible postoperative complications.
- about the material composition of the implant.
- that implants may respond to metal detectors during security checks (e.g. at airports) and carrying an implant ID as a proof is recommended.
- that implants may interact with medical imaging technique (e. g. MRI)

22. Circumstances that can interfere with the Success of an Operation

- Severe osteoporosis
- Severe deformities
- Local bone tumors
- Systemic diseases
- Metabolic disorders
- Case history of infections and falls
- Drug dependency or abuse, including excessive alcohol and nicotine consumption
- Obesity
- Mental disorders or neuromuscular diseases
- Heavy physical activities associated with strong vibrations
- Hypersensitivities

23. Warnings / Precautions

- The reuse of LINK single-use products is not permissible.
- Implants must be handled with great care and should not be modified or changed, even the smallest scratches and damages can considerably impair their stability or performance. Damaged implants are not permitted to be used.
- Surfaces provided for the connection of modular prosthetic components (cone, pins, screws) must not be damaged and may need to be cleaned with sterile liquid and dried before being joined together, so that neither blood nor any other coating impairs any of the connections, which could compromise the reliability of the connection.
- Do not manipulate or misuse instruments. We do not accept liability for products that have been modified, subjected to unintended use, or used improperly.
- For the processing of LINK instruments, it is presumed that the personnel have technical knowledge level I (Germany) and in other countries technical knowledge and expertise.
- Medical devices that are sent in for servicing must be processed beforehand in such a way that they cannot constitute a hazard to third parties.
- Products made of plastic (e.g. polyamide (PA), polyethylene (PE), polyoxymethylene (POM), ultra-high molecular weight polyethylene (UHMWPE)) may not be localizable using external imaging procedures.

24. Preoperative Planning

Preoperative planning provides important information to identify the appropriate implant system and select the components of a system. Make sure that all components required for the operation are laid out and ready in the operating room. Trial implants to verify proper fit (where applicable) and additional implants should be kept at the ready, in case other sizes are needed or the intended implant cannot be used. All LINK instruments necessary for the implantation must be on-hand, sterilized and intact.

If prosthesis implantation is indicated, then it must be taken into consideration, along with the overall circumstances of the patient:

- that all non-surgical and surgical treatment alternatives for the joint disease have been considered
- that artificial joint replacement performance is categorically inferior to natural joint performance, and an indication-related improvement in the preoperative condition is the only aim here
- that proper selection, placement and fixation of the devices are decisive factors, which will determine the life of the implant.
- that an artificial joint may loosen due to stress, wear and tear, and infection, or luxation or dislocation may occur
- that revision surgery, which under certain circumstances may exclude the possibility of restoring joint function, may be necessary due to loosening of the implant
- that the patient consents to undergo the operation and accepts the risks involved
- that if load-transferring bone structures are damaged, then the loosening of the components, bone and implant fractures, as well as other serious complications cannot be ruled out
- that if the patient is suspected of having allergies and tests positive on the applicable tests, then the patient's foreign body sensitivities (material tolerances) must be examined
- that acute and chronic infections, local and systemic, may compromise the successful implantation, so pre-operative microbiological analysis is recommended

Generally, the mechanical failure or fracture of an implant is a rare exception. However, this cannot be excluded with absolute certainty despite the sound structure of the implant.

This may be due to stress on the implant and prosthesis as the result of a fall or accident, among other things.

If the bone area where the implant is anchored is altered in such a way that the prosthesis is no longer able to withstand normal stress and an area of the prosthesis becomes subject to a stress imbalance, then a mechanical failure of the implant system may result. Such stress imbalances may also occur if the anchoring elements of the implants are required to form a bridge over larger bone deficiencies without optimal reinforcement of the bone. It is recommended that the implant with the largest possible anchoring elements be used. Proper preparation for surgical procedures also includes the functional testing of implants and instruments prior to use.

For definitive identification information on the product such as system compatibility, article number, material and shelf life, refer to the identification on the implant and / or the packaging. You should also take advantage of the training courses and printed materials provided for your information. To learn more, please contact the Waldemar Link GmbH & Co. KG sales office or your field representative.

25. Handling

All implant components are supplied sterile as single-use devices in individual packages. The implant components are sterilized by gamma sterilization, at least 25 kGy.

By contrast, components made of highly cross-linked polyethylene or highly cross-linked polyethylene with vitamin E are sterilized with ethylene oxide (ETO).

Implants should always be stored in their unopened protective packaging. Examine the packaging for damage before using the implant. Damaged packaging can have an adverse effect on both the sterility of the device as well as the proper performance of the implant, such that the device may no longer be used.

- Check the use by date on the implants. Implants with expired use by dates are no longer permitted to be used for implantation!
- After opening the package, check that model and size of the implant are matching with the information printed on the package labeling.
- Observe the pertinent standards for the aseptic handling of devices during and after removal of the implant from the packaging.
- When removing the packaging, make a record of the batch or serial numbers on the label, since this information is decisive for batch tracing. Self-adhesive labels with this information are enclosed with every package for your convenience.

Instruments must always be treated with care, this particularly applies during transport, cleaning, maintenance, sterilization, and storage. The sterile status of the instruments depends, inter alia, on the sterile items packaging and the prevailing storage conditions and must be established together with the operator's hygiene officer on a case-by-case basis. Direct sunlight must be avoided. Improper handling and care, as well as unintended use, can lead to premature wear or damage.

Devices made of plastics (e.g. PP-H) may not be located by means of an external imaging device.

26. Intraoperative Use

Please refer to the relevant surgical technique associated with the system for information about the intraoperative use of the system.

27. Postoperative

In addition to movement and muscle training, special attention must be paid to carefully instructing the patient during the postoperative phase.

Physician-supervised postoperative monitoring of healing progress is recommended. Where applicable, patients should also be advised on how to avoid overstraining themselves.

Follow-up examinations should be carried out regularly or immediately if symptoms occur.

28. Notes on MRI and CT Examination Procedures

Our implants were not evaluated for safety and compatibility for MRI and CT examination procedures.

In the case of our metallic implants and implant components, MRI examinations pose potential risks to the patient due to possible heating and migration of the implants or implant components.

Likewise, there is a potential risk of artifact formation in MRI and CT examinations of our metallic implants and implant components.

The probability of occurrence and the extent of the potential risks mentioned depend on the type of device used, its device parameters and the sequences used.

Always follow the instructions in the operating instructions of the manufacturer of the device used for the imaging.

The selection of the imaging examination procedure and the assessment of possible side effects is the responsibility of the examining physician.

The examining physician must take into account the individual condition of the patient and other diagnostic methods.

29. Explantation of Implants / Revision Surgery

Please refer to the relevant surgical technique associated with the system for information about explantation of the implants and revision surgery.

30. Disposal

Packaging and system components to be discarded must be handled in compliance with your national and local regulations for hospital disposal.

31. Instruments

Please refer to the description in the reprocessing instructions H50 for the:

- initial use
- performance test
- maintenance
- manual cleaning
- cleaning in a washer-disinfector
- reprocessing
- sterilization
- servicing
- transport

32. Link to Summary of Safety and Clinical Performance

The Summary of Safety and Clinical Performance of this product can be found on the EUDAMED database.

<https://ec.europa.eu/tools/eudamed/>

Please, search by UDI-DI number of this product. The UDI-DI number is shown on the product label.

33. Requests

Requests of any kind should be directed to Waldemar Link GmbH & Co. KG (see contact information in this document).

34. Complaints about our Products

All complaints must be addressed to Waldemar Link GmbH & Co. KG at: complaint@link-ortho.com

In the event of a complaint, the name or reference number of the corresponding component should be specified with the serial number (SN) or the lot number (LOT), your name, and your contact address. The reason for the complaint should be given in brief.

35. Report of serious incidents

Any serious incident that occurs in relation to the device must be reported to the manufacturer and the authority responsible for your location.

Waldemar Link GmbH & Co. KG and/or other corporate affiliated entities own, use or have applied for the following trademarks in many jurisdictions: LINK, BiMobile, SP II, Modell Lubinus, E-Dur, EndoDur, T.O.P. II, BetaCup, CombiCup PF, CombiCup SC, CombiCup R, MobileLink, C.F.P., LCU, SP-CL, LCP, MIT-H, Endo-Model, Endo-Model SL, MP, MEGASYSTEM-C, GEMINI SL, SPAR-K, LCK, Link OptiStem, HX, TiCaP, X-LINKed, PorAg, LINK PorEx, BiPorEx, PorEx-Z, TrabecuLink, Tilastan, customLINK, RescueSleeve, Stactip, VACUCAST.

Other trademarks and trade names may be used in this document to refer to either the entities claiming the marks and/or names or their products and are the property of their respective owners.

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