







## TrabecuLink Tibial Cones

Stable – Elastic – Versatile

**CE0482**

Explanation of Pictograms			
	Manufacturer		Article number
	Material (number)		Product meets the applicable requirements, which are regulated in the EU harmonization legislation for the affixing of the CE marking.

The dynamic **LINK Tibial Cones** are an attractive solution for cementless restoration of bone defects<sup>10</sup> and to provide additional support for the prosthesis if there is bone loss in the proximal tibia. The combination of the **dynamic design**<sup>5,6</sup> of the cones and the biocompatible material **Tilastan-E**<sup>11,12</sup> is ideal for ensuring stable, long-lasting fixation and successful bone regeneration.

The **3-dimensional TrabecuLink structure**, with its pore size, porosity and structure depth, also provides an excellent basis for promoting osteoconduction and microvascularization, taking into account the requirements for the structure-covering protein layer (fibronectin - vitronectin - fibrinogen).<sup>1,2</sup>

Tibial Cones can be used in combination with the long-established **LINK Endo-Model knee family** in a wide range of sizes and versions. The choice of sizes corresponds to the dimensions of the hinged knee prostheses.


- **Stable** – with cementless fixation
- **Elastic** – due to integral bending axes in the inner metal wall
- **Versatile** – for a broad range of solutions<sup>7</sup>






**4 Sizes and 4 Versions**


- XS, S, M, L
- Full, right-half, left-half, half




**Cementless Implantation  
in the Bone**  
for bone regeneration




**Extraction Pockets**  
for gripping with instruments



**Inner Metal Wall**  
with parallel 3 mm notches



**Spring Effect**  
provided by elastic bending axes  
and the compensator (for ML and  
AP elasticity)



**TrabecuLink Structure**

- Pore size: 610-820 µm
- Structure depth: 2 mm
- Porosity: 70%
- Bone reaction-friendly

1 Real TrabecuLink® structure  
2 Geometric TrabecuLink® structure

## Stable – in metaphyseal fixation<sup>9,13</sup>

- Reinforcement of the bone structure in cases of tibial, metaphyseal bone defects
- High primary stability in the metaphyseal region, both for the cone itself and for the tibial component cemented in the cone
- Cementless interface with the bone for bone regeneration

## Elastic – due to integral bending axes in the inner metal wall

- Mechanical compression promotes bone regeneration<sup>5,6</sup>
- Bending axes for adaptation to bone surfaces: oriented vertically in the cone
- Good fit ensured by structural elasticity, which also facilitates insertion of the tibial cones
- Spring effect for easier intraoperative positioning

## Versatile – for a broad range of solutions<sup>7</sup>

- Can be combined with all the tibial components of the LINK Endo-Model knee family
- Sizes correspond to the sizes of the hinged knee prostheses
- Customized models can be manufactured

## Protective – due to inner metal wall

- Prevents penetration of bone cement into the TrabecuLink structure
- Reliable cement fixation by means of specially positioned “notches” (revision-friendly)

## Environmentally friendly<sup>3,8</sup>

- Resource-saving manufacturing of proven Titanium alloy






Example of a patient-specific custom-made implant



## TrabecuLink

3-dimensional structure – for optimal bone ongrowth

- Pore geometry (porosity: 70%, pore size: 610-820  $\mu\text{m}$ , structure depth: 2 mm) ensures excellent cell ongrowth<sup>1,2,4</sup>



## Biomaterials

for joint prostheses

Biomaterials for joint prostheses should not be cytotoxic and, if intended for cementless anchorage, they should promote the deposition of bone precursor cells. On this basis, Steinemann<sup>4</sup> conducted tests with samples of widely differing materials in cultures with fibroblasts and osteoblasts. The result shows that titanium and zirconium were the only implant materials that supported osteoblast proliferation. Niobium and tantalum, in contrast, greatly hindered growth.<sup>4</sup>

Further information on this can be found in the special issue **DirectLINK 1-2017 “Periprosthetic Infections”** at <https://www.linkorthopaedics.com>.



1 Sizing the implant



2 Insertion of the trial cone



- 3 Trial reduction:
- Trial tibial component
  - Trial stems
  - Trial washers as required
  - Trial cone
  - Trial femoral component



4 Insertion of the final TrabecuLink Tibial Cone



5 TrabecuLink Tibial Cone in situ




6 Cement application and insertion of the tibial component

**Note:** The surgical procedures for Endo-Model implants are described in the implant-specific surgical techniques.



## TrabecuLink Tibial Cones

For use with Endo-Model tibial components standard and modular version,  
Material: Tilastan – E (TiAl6V4)



### Combination Chart:

Sizes of Knee Prostheses

	XS	S	M	L
XS				
S				
M				
L				



Combination for the entire Endo-Model knee family (Endo-Model, Endo-Model –M, Endo-Model SL\*)



Combination only with Endo-Model and Endo-Model –M




Only after prior use of all trial components and 5 mm deeper positioning (Endo-Model, Endo-Model –M)



This combination is not permitted

\* The geometry of the Endo-Model SL tibial components only allows a combination with tibial cones in the same size.



**References (general)**

- <sup>1</sup> Cecile M. Bidan, Krishna P. Kommareddy, Monika Rumpler, Philip Kollmannsberger, Yves J.M. Brechet, Peter Fratzl, John W.C. Dunlop. et al.;  
**How Linear Tension Converts to Curvature: Geometric Control of Bone Tissue Growth;**  
*PLoS ONE* 7(5): e36336. <https://doi.org/10.1371/journal.pone.0036336> (2012)
- <sup>2</sup> Pascal Joly, Georg N. Duda, Martin Schöne, Petra B. Welzel, Uwe Freudenberg, Carsten Werner, Ansgar Petersen, et al.;  
**Geometry-Driven Cell Organization Determines Tissue Growth in Scaffold Pores: Consequences for Fibronectin Organization;**  
*PLoS ONE* 8(9): e73545. <https://doi.org/10.1371/journal.pone.0073545> (2013)
- <sup>3</sup> Dr. Malte Drobe, Franziska Killiches;  
**Vorkommen und Produktion mineralischer Rohstoffe – ein Ländervergleich;**  
*Bundesanstalt für Geowissenschaften und Rohstoffe Hannover;*  
[http://www.bgr.bund.de/DE/Themen/Min\\_rohstoffe/Downloads/studie\\_rohstoffwirtschaftliche\\_einordnung\\_2014.pdf?\\_\\_blob=publicationFile&v=4](http://www.bgr.bund.de/DE/Themen/Min_rohstoffe/Downloads/studie_rohstoffwirtschaftliche_einordnung_2014.pdf?__blob=publicationFile&v=4) (2014)
- <sup>4</sup> Steinemann SG;  
**Compatibility of Titanium in Soft and Hard Tissue – The Ultimate is Osseointegration;**  
*Materials for Medical Engineering, WILEY-VCH, Volume 2, Page 199-203*
- <sup>5</sup> Gerald Küntscher; **Praxis der Marknagelung;** *Friedrich-Karl Schattauer-Verlag (1962)*
- <sup>6</sup> R. Texhammer, C. Colton et al.; **AO-Instrumente und Implantate (Technisches Handbuch);**  
*Springer Verlag, 2. Auflage, S.25 (2011)*
- <sup>7</sup> Gabriele Panegrossi, corresponding author Marco Ceretti, Matteo Papalia, Filippo Casella, Fabio Favetti, and Francesco Falez;  
**Bone Loss Management in Total Knee Revision Surgery;** *Int Orthop.* 2014 Feb; 38(2): 419–427;  
<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC3923937/> (2014)
- <sup>8</sup> **Conflict Minerals: MEPs Secure Mandatory Due Diligence for Importers;**  
*Press release - External/international trade – 22-11-2016 - 19:07;* <http://www.europarl.europa.eu/news/en/news-room/20161122IPR52536/conflict-minerals-meps-secure-mandatory-due-diligence-for-importers> (2016)
- <sup>9</sup> Henricson A, Linder L, Nilsson KG.;  
**A Trabecular Metal Tibial Component in Total Knee Replacement in Patients Younger than 60 Years: a Two-year Radiostereophotogrammetric Analysis;** *J Bone Joint Surg Br.* 2008;90:1585–1593. doi: 10.1302/0301-620X.90B12.20797 (2008)
- <sup>10</sup> P. K. Sculco, M. P. Abdel, A. D. Hanssen, D. G. Lewallen;  
**The Management of Bone Loss in Revision Total Knee Arthroplasty;**  
*Bone Joint J* 2016;98-B(1 Suppl A):120–4 (2016)
- <sup>11</sup> Peter Heintl, Lenka Müller, Carolin Körner, Robert F. Singera, Frank A. Müllerb;  
**Cellular Ti–6Al–4V Structures with interconnected Macro Porosity for Bone Implants Fabricated by Selective Electron Beam Melting;**  
*Acta Biomaterialia* Volume 4, Issue 5, September 2008, Pages 1536–1544 (2008)
- <sup>12</sup> Hong Wang, Bingjing Zhao, Changkui Liu, Chao Wang, Xinying Tan, Min Hu;  
**A Comparison of Biocompatibility of a Titanium Alloy Fabricated by Electron Beam;**  
*PLOS ONE* | DOI:10.1371/journal.pone.0158513 July 8 2016, (2016)
- <sup>13</sup> Ivan De Martino, Vincenzo De Santis, Peter K Sculco, Rocco D’Apolito, Joseph B Assini, Giorgio Gasparini;  
**Tantalum Cones Provide Durable Mid-Term Fixation in Revision TKA;**  
*Clin Orthop Relat Res* 473 (10), 3176-3182 (2015)

Please note the following regarding the use of our implants:

**1. Choosing the right implant is very important.**

The size and shape of the human bone determines the size and shape of the implant and also limits the load capacity. Implants are not designed to withstand unlimited physical stress. Demands should not exceed normal functional loads.

**2. Correct handling of the implant is very important.**

Under no circumstances should the shape of a finished implant be altered, as this shortens its life span. Our implants must not be combined with implants from other manufacturers. The instruments indicated in the Surgical Technique must be used to ensure safe implantation of the components.

**3. Implants must not be reused.**

Implants are supplied sterile and are intended for single use only. Used implants must not be used again.

**4. After-treatment is also very important.**

The patient must be informed of the limitations of the implant. The load capacity of an implant cannot compare with that of healthy bone!

**5. Unless otherwise indicated, implants are supplied in sterile packaging.**

Note the following conditions for storage of packaged implants:

- Avoid extreme or sudden changes in temperature.
- Sterile implants in their original, intact protective packaging may be stored in permanent buildings up until the "Use by" date indicated on the packaging.
- They must not be exposed to frost, dampness or direct sunlight, or mechanical damage.
- Implants may be stored in their original packaging for up to 5 years after the date of manufacture. The "Use by" date is indicated on the product label.
- Do not use an implant if the packaging is damaged.

**6. Traceability is important.**

Please use the documentation stickers provided to ensure traceability.

**7. Further information** on the material composition is available on request from the manufacturer.

**Follow the instructions for use!**

### Waldemar Link GmbH & Co. KG, Hamburg


All content in this catalog, including text, pictures and data, is protected by law. Every instance of use, whether in part or in whole and which is not permitted by law, is subject to our prior consent. In particular, this applies to the reproduction, editing, translation, publishing, saving, processing, or passing on of content stored in databases or other electronic media and systems, in any manner or form. The information in the catalogs is solely intended to describe the products and does not constitute a guarantee.

The Surgical Technique described has been written to the best of our knowledge and belief, but it does not relieve the surgeon of his/her responsibility to duly consider the particularities of each individual case.

Products shown in this document may not be available in your country. The product availability is subject to the approval and/or registration regulations of the respective country. Please contact Waldemar Link GmbH & Co. KG if you have questions about the availability of LINK products in your country.

Waldemar Link GmbH & Co. KG and/or other corporate affiliated entities own, use or have applied for the following trademarks in many jurisdictions: LINK, BiMobile, SP II, Modell Lubinus, E-Dur, EndoDur, T.O.P. II, BetaCup, CombiCup PF, CombiCup SC, CombiCup R, MobileLink, C.F.P., LCU, SP-CL, LCP, MIT-H, Endo-Model, Endo-Model SL, MP, MEGASYSTEM-C, GEMINI SL, SPAR-K, LCK, Link OptiStem, HX, TiCaP, X-LINKed, PorAg, LINK PorEx, BiPorEx, PorEx-Z, TrabecuLink, Tilastan, customLINK, RescueSleeve, Stactip, VACUCAST.

Other trademarks and trade names may be used in this document to refer to either the entities claiming the marks and/or names or their products and are the property of their respective owners.

 Waldemar Link GmbH & Co. KG

Barkhausenweg 10 · 22339 Hamburg · Germany

Phone +49 40 53995-0 · [info@linkhh.de](mailto:info@linkhh.de)

[www.linkorthopaedics.com](http://www.linkorthopaedics.com)

**LINK**<sup>®</sup>

